

Klerwipe™ Nonwoven Dry Wipe

PHARMACEUTICAL

General purpose cleanroom wipe. Made from a binder free 55% cellulose/45% polyester hydroentangled nonwoven blend. Recommended for use in ISO Class 5 and above environments.

The validated sterile version is ideal for cleaning applications during aseptic processing activity. Excellent for use in Primary Engineering Controls (PEC) such as laminar air hoods, biological safety cabinets and compounding aseptic isolators where sterility of end product must be assured.

Nonwoven

FEATURES

- Autoclavable
- Double bagged
- High absorbency and durable
- Chemically compatible with disinfectants and cleaning solutions
- Nonwoven dry wipers are more absorbent than knitted wipes. Binder free means no residue are left on critical surfaces
- An excellent value combined with high performance cleaning characteristics

Validated Sterile

FEATURES

- Gamma irradiated to a 10⁻⁶ Sterility Assurance Level (SAL)
- Certificates of Sterility available
 - Expiration date
 - Radiation dosage
 - Lot information
 - Meets cGMP requirements for traceability
- Double bagged
- Header with Package Integrity testing

APPLICATIONS

- General wiping
- Component preparation
- Wash areas
- Sterile gowning rooms
- Applying and removing disinfecting liquids



WORLDWIDE HEADQUARTERS

1 Ecolab Place
St. Paul, MN 55102
USA

www.ecolab.com/lifesciences
www.ecolab.com/pharmaceutical

© 2019 Ecolab USA Inc. All rights reserved.

ORDER INFORMATION

NONWOVEN

PART #	53000344
SIZE	9X9"
WIPES/PACK	300
PACKS/CASE	20

PART #	53000346
SIZE	12X12"
WIPES/PACK	150
PACKS/CASE	20

VALIDATED STERILE

PART #	53000342
SIZE	9X9"
WIPES/PACK	20
PACKS/CASE	80

PART #	53000343
SIZE	9X9"
WIPES/PACK	50
PACKS/CASE	20

PART #	53000345
SIZE	12X12"
WIPES/PACK	20
PACKS/CASE	25

FOR ADDITIONAL INFORMATION, PLEASE CONTACT YOUR ACCOUNT MANAGER OR VISIT OUR WEBSITE

ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

ECOLAB®
Everywhere It Matters.™